# **Electricité**

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# ***Electricity***

*Electromagnetism*

**Ref :  
292 060**

**Rails de Laplace**


***Laplace's rails***

Français – p 1

English – p 3

Version : 9003

## 1 Description

### 1.1 Composition

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
## 1 Description

### 1.1 Composition

- 1) 2 aluminium rails fixed on a supporting base and connected to 2 safety banana plugs
- 2) 1 movable bar fitted with a marked disk
- 3) 2 red and blue markers to indicate the poles of the magnet and direction of the magnetic field

### 1.2 Additional equipment required (not supplied)

- 1 U-shaped magnet **Part no. 263 006**
- R30 30V-5A adjustable power supply **Part no. 281 090**

#### **IMPORTANT**


**If using the device for the first time or if the device has not been used for a long period, it may be necessary to clean the aluminium parts with a fine sandpaper. Otherwise, there may be a bad contact between the rails and the movable bar may not move or may move very little.**

## 2 Setting up

Place the device horizontally on a table. Place the U-shaped magnet between the rails and install the blue and red markers to indicate the direction and line of the magnetic field to the class.

Place the movable bar on the rails, perpendicular to these and the marked disk in such a way that the movement is visible to the entire class.

Connect the banana plugs to the power supply and pass a current of at least 2.5 A in the circuit. The movable rod then is subjected to Laplace's force that makes it roll on the rails as shown in the diagram below.

### 3 Technical characteristics

- **Current intensity for starting the movement:** 2.5 A
- **Rail and movable bar:** aluminium
- **Connections:** safety banana sockets  $\varnothing$  4 mm

### 4 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

## Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

### Free service \*

**+ 33 (0)2 32 29 40 50\*\***

\* Call cost not included

\*\* Only for call from foreign countries